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Warning: Please read this instruction manual carefully before using the instrument and keep it for future use.

A. IMPORTANT INFORMATION

1. DO NOT use it to or around the mucosa, wounds or eyes. Children MUST use it under the guidance of adults.
2. Infants, people who have no consciousness or sense, and people who are insensitive to temperature are FORBIDDEN to use this product.
3. This product has temperature autocontrol function. There is a risk of doing harm to your skin if it is applied to healthy skin or when the itching feeling has disappeared.

This device uses infrared heat to quickly and durably relieve the pain and itching that follows an insect bite.

The concentrated heat produced by the infrared chip naturally breaks down the proteins and thus eliminates the toxins causing the discomfort associated with the sting.

To ensure user safety, this product is chemical free, additive free, and uses a medical grade thermal conductivity material.

B. SAFETY INSTRUCTIONS

1. The device is intended exclusively for private and not for commercial use. Only use this device as describe in this instruction manual.
2. Never try to repair it by yourself. It must be repaired by the manufacturer or its service agent or a similarly qualified person in order to avoid an electric shock.
3. Keep it away from heat sources such as radiators, to avoid deformation of plastic parts.

C. HOW TO USE IT

1) Insert or replace batteries

- Remove the battery cover on the back of the unit.
- Insert 2 AA batteries observing the correct polarity indicated.
- Replace the cover.

The red LED's lighting on means that the batteries run out and need to be replaced.

2) Operation

- Let the head of physical therapy close to the itching or bite position (Must be sure to turn on the power after applying the physical therapy cover close to the bite).
- Press the button to turn on the device, with the flash of blue LED illustrating the device in normal operation.
- The device can be switched from grade 1 to 5 to get the different appropriate temperature (Generally children fit grade 1 to 2, while adults fit grade 2 to 5. Different people and each body part have different feelings of temperature so the heat grade should depend on the users' actual situation).
- Keep using 2~3 seconds after the itching gradually disappearing and then press the button again to turn off the device.

D. CARE AND MAINTENANCE

Clean the product only with a dry or slightly damp, lint-free cloth.
Do not use abrasive cleaners to clean the unit.

E. FEATURES

Material : ABS

Power by 2xAA batteries (not included)

Heat: 47/49°C

The European directive 2012/19/EC about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capacities, or by persons lacking in experience or knowledge, unless they are supervised or instructed on the use of this apparatus by a person responsible for their safety. Children should be supervised to ensure they do not play with this appliance.

Precision and precautions about the battery use

- Throwing away a used battery in the environment or with household waste pollutes and avoids the recuperation of reusable materials.
- It is important to limit the battery consumption and to respect the following rules:
 - Favor alkaline battery (last longer than saline battery) and when it's possible use rechargeable battery.
 - Leave used batteries and accumulators in adapted containers that you can find in some retailers. In that way, metals don't pollute the environment and the health.
- Battery has to be inserted by respecting the polarity show on the apparel and the battery. An incorrect positioning can damage the apparel, cause leaking and even provoke fire and the explosion of the battery.
- To assure a good functioning, battery has to be in a good condition. If there is a problem, put new battery.
- Never try to recharge battery not made for this. It could leak, warm up, cause fire or explode.
- Replace all the battery at the same time. Never mix up alkaline battery with saline battery.
- Used battery has to be removed.
- Also, remove the battery if you're not using the apparel for a long time.
- Never try to short-circuit the terminal of a battery.
- Never throw away the battery in fire because they can explode.
- Take out the accumulator from the apparel before reload it.
- We recommend that an adult watch if a child change the battery to respect those rules.
- If you swallow a battery, go immediately to your doctor or the hospital.
- For the best ongoing performance, we recommend using only fresh, high-quality, alkaline batteries.

Importé par / Importiert von / Geïmporteerd door / Imported by :
CAPTELEC, 59170 Croix - FRANCE

Carton et emballage papier à trier / Kartons und papierverpackung recyceln / Verpakkingsmateriaal in karton en papier te sorteren / Recycle carton and paper packaging